
ECO pH Bedienungsanleitung

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
Front panel and display

Technical data

Range	0 to 14 pH, -5 to 60°C (23 to 140°F)
Resolution	0,01 pH, 0,1°C (0,1°F)
Accuracy	± 0,01 pH, ± 0,5°C (± 1°F)
Calibration	automatic 1- or 2-point calibration
Temperature compensation	automatic
Battery	4 x 1,5V button type, approx. 100 hours of use
Auto-off	after 8 minutes of non-use
Envirment	-5 to 50°C (23 to 122°F), RH 100%
Housing	waterproof
Dimension	175,5 x 39 x 23 mm
Weight	100 g
Warranty	1 year, 3 month electrode

Operation

Turn on and measurement

- Press and hold **ON/OFF** until the display lights up.
- Remove the protective cap of the electrode. The tester is now ready for operation.
- Put the electrode into the sample. The stability indicator shows a stable measurement result.
- Freeze the display with **SET/HOLD**. Press any button to return to measurement mode.
- To turn off the tester press **ON/OFF** until the display shows **OFF**.
- Clean the electrode after each measurement session with distilled water. Put a few drops calibration solution pH 4,01 in the protective cap and put on the cap.

*Do **never** store the electrode in distilled water.*

Calibration

Calibrate the ECO pH-Tester, if

- the electrode has been replaced
- min. once per month
- to ensure accurate results
- after measuring aggressive chemicals

- Press **ON/OFF** to turn the tester on.
- After a few seconds press and hold **ON/OFF**, until the display shows **OFF** followed by **CAL**.
- When **CAL** is shown release **ON/OFF**.
- The display shows **pH 7.01 USE**.
- Put the electrode in the calibration solution pH 7.01. The display shows **OK1**.
- To finish the **1-point calibration**, press **ON/OFF**, the tester returns to the measurement mode.
 For a **2-point-calibration**, clean the electrode and put it in the calibration solution pH 4.01 or 10.01. The display shows **pH 4.01 USE**.
- After about 12 seconds the display shows **OK 2**, when the calibration is successful. The tester returns to the measurement mode.

Maintenance

Change the electrode

- Unscrew the retaining ring and put it off.
- Pull off the old electrode body.
- Put on the new electrode body.
- Screw the retaining ring on.
- Calibrate the electrode.

Change the battery

- When the display shows the battery icon, change the batteries.
- Unscrew the retaining ring and put it off.
- Pull off the electrode body.
- Pull out the battery compartment gently.
- Change the four 1,5V button batteries, pay attention to the polarity.
- Press the battery compartment carefully back into the ECO pH-Tester.