Datasheet HCA-S-400M-IN

400 MHz Photoreceiver with InGaAs-PIN Photodiode



The picture shows model HCA-S-400M-IN-FS. The photoreceiver will be delivered without post holder and post.


Features

- InGaAs-PIN photodiode
- Bandwidth DC 400 MHz
- Amplifier transimpedance gain 5.0 × 103 V/A
- Max. conversion gain 4.8×10^3 V/W @ 1550 nm
- Spectral range 900 1700 nm
- Free-space input 1.035"-40 threaded, alternatively 25 mm diameter unthreaded
- Fiber optic input available as permanently mounted FC-input
- UNC 8-32 and M4 tapped holes for mounting on standard posts with metric and imperial thread

Applications

- Spectroscopy
- · Fast pulse and transient measurements
- Optical triggering
- Optical front-end for oscilloscopes, A/D converters and HF lock-in amplifiers

Block Diagram



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E 🍑 T O

HCA-S-400M-IN_R6/TH,JMa/04MAY2022 Page 1 of 10

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Available Versions

HCA-S-400M-IN-FST



Picture shows 1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm) 1.035"-40 threaded flange for free space applications. Compatible with many optical standard accessories.

HCA-S-400M-IN-FS



Picture shows unthreaded flange with 25 mm diameter

25 mm dia. unthreaded flange for free space applications. Compatible with many optical standard accessories.

HCA-S-400M-IN-FC



Fix/permanent FC fiber connector for high coupling efficiency and excellent conversion gain accuracy.

Related Models

HCA-S-400M-SI-FST

Si-PIN. Ø 0.8 mm. 320 − 1000 nm

free space input, 1.035"-40 threaded flange

HCA-S-400M-SI-FS

 $Si-PIN,\emptyset$ 0.8 mm, 320-1000 nm

free space input, 25 mm dia. unthreaded flange

Si-PIN,∅ 0.8 mm, 320 – 1000 nm

HCA-S-400M-SI-FC

FC fiber connector (fix/permanent)

Si-PIN, Ø 0.8 mm, 320 − 1000 nm

HCA-S-400M-SI-SMA

SMA fiber connector (fix/permanent)

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

Datasheet

HCA-S-400M-IN

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Available Accessories

PRA-PAP



Alternative mounting option: Post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, PWPR, HCA-S and LCA-S.

PS-15-25-L



Power supply Input: 100 - 240 VAC Output: ±15 VDC

Specifications

Test conditions

 $V_S = \pm 15 \text{ V}$, $T_A = 25 \,^{\circ}\text{C}$, output load impedance 50 Ω , warm-up 20 minutes (min. 10 minutes recommended)

Gain

Transimpedance gain Gain accuracy Conversion gain

 5.0×10^3 V/A (@ output load 50 Ω)

±1 % (electrical) 4.8×10^{3} V/W typ. (@ 1550 nm, output load 50 Ω)

Frequency Response

Lower cut-off frequency Upper cut-off frequency (-3 dB)

400 MHz (±15 %)

Gain flatness

 $\pm 1 dB$

Time Response

Rise/fall time (10 % - 90 %)

1.0 ns

Noise equivalent power (NEP) Input

Optical saturation power Input offset compensation range

24 pW/√Hz (@ 1550 nm, 100 MHz) 200 uW (for linear amplification, @ 1550 nm) ±200 μA, adjustable by offset potentiometer

Detector

Detector Active area (FS/FST version) Active area (FC version)

InGaAs-PIN photodiode Ø 0.3 mm integrated ball lens

suitable for fibers up to 62.5 µm core diameter

Spectral range Max. sensitivity

900 - 1700 nm 0.95 A/W typ. (@ 1550 nm)

Output

Output voltage range

 $\pm 1.0 \text{ V } (@ 50 \Omega \text{ output load})$

for linear operation and low harmonic distortion $\pm 1.5 \text{ V}$ (@ 50 Ω output load)

Max. output voltage range Output impedance Output noise

50 Ω (terminate with 50 Ω load)

3 mV_{RMS} (20 mV_{PP}) typ. (@ 50 Ω load, no signal on detector, measurement bandwidth 1.5 GHz)

Input Flange

Material

1.4305 stainless steel, nickel-plated (FST flange) AlMg4.5Mn, nickel-plated (FS flange)

Coupler Ring

(FST version only)

Material

1.4305 stainless steel, glass bead blasted

Power Supply

Supply voltage Supply current ±15 V (±14.5 V ... ±16.5 V)

±55 mA (depends on operating conditions,

recommended power supply capability min. ±150 mA)

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

0

Datasheet


HCA-S-400M-IN

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Specifications (continued)			
Case	Weight	209 g (0.46 lbs) HCA-S-400M-IN-FST incl. coupler ring 196 g (0.43 lbs) HCA-S-400M-IN-FS 188 g (0.41 lbs) HCA-S-400M-IN-FC AIMg4.5Mn, nickel-plated	
	Material		
Temperature Range	Storage temperature Operating temperature	-30 °C +85 °C 0 °C +60 °C	
Absolute Maximum Ratings	Optical input power (CW) Power supply voltage	10 mW ±20 V	
Connectors	Input	HCA-S-400M-IN-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories
		HCA-S-400M-IN-FS	25 mm dia. unthreaded flange for free space applications
		HCA-S-400M-IN-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible)
	Output	BNC jack (female)	
	Power supply LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)		
		PIN 2 O O O	PIN 1 +Vs Pin 1: +15 V Pin 2: -15 V Pin 3: GND
Scope of Delivery	HCA-S-400M-IN, internally threaded coupler ring (FST version only), LEMO® 3-pin connector, datasheet, transport package		
Ordering Information	HCA-S-400M-IN-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories.	
	HCA-S-400M-IN-FS	25 mm dia. unthreaded flange for free space applications.	
	HCA-S-400M-IN-FC	FC fiber optic connector (fix/permanent, FC/PC and FC/APC compatible).	

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E MT O




PD-HCA-5-400M-IN-DW_KU I

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

E M T O

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Typical Performance Characteristics (continued)




PD-HCA-S-400M-IN-noise R0

Note: spectral noise data is measured at the amplifier output with no signal on the photodiode. To determine the spectral input noise divide the measured output noise by the amplifier conversion gain.

Conversion gain (V/W) = amplifier gain $(V/A) \times$ photo sensitivity (A/W).

Marker frequency output noise resulting input noise (NEP) $\frac{1}{100 \text{ MHz}} \frac{112 \text{ nV/}\sqrt{\text{Hz}}}{12 \text{ pV/}\sqrt{\text{Hz}}} \frac{24 \text{ pW/}\sqrt{\text{Hz}}}{1550 \text{ nm}}$

Pulse response to square wave input signal (with 16 times averaging)




PD-HCA-S-400M-IN-pulse-2ns_R01

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E 🍑 T O


400 MHz Photoreceiver with InGaAs-PIN Photodiode

Typical Performance Characteristics (continued) Large signal response output signal for 100 MHz, 210 μW modulated optical input signal (with 4 times averaging)



PD-HCA-S-400M-IN-large-sinus_R01

Small signal response output signal for 3 µW modulated optical input signal, 1 MHz square wave, without averaging




SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

0


400 MHz Photoreceiver with InGaAs-PIN Photodiode

Typical Performance Characteristics (continued) Small signal response output signal for 3 μ W modulated optical input signal, 1 MHz square wave, with 64 times averaging



Dimensions

HCA-S-400M-IN-FST (1.035"-40 threaded free space input)



SOPHISTICATED TOOLS FOR SIGNAL RECOVERY


F E M T O


Page 8 of 10

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Dimensions (continued)

HCA-S-400M-IN-FS (25 mm dia. unthreaded free space input)





DZ-HCA-S-X00-IN_FS_R1

all dimensions in mm unless otherwise noted

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O


Datasheet


HCA-S-400M-IN

400 MHz Photoreceiver with InGaAs-PIN Photodiode

Dimensions (continued)

HCA-S-400M-IN-FC (FC fiber optic connector)





DZ-HCA-S-XX-XX_FC_R1

all dimensions in mm unless otherwise noted

FEMTO Messtechnik GmbH Klosterstr. 64 10179 Berlin · Germany Phone: +49 30 280 4711-0 Fax: +49 30 280 4711-11 Email: info@femto.de www.femto.de

Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany

SOPHISTICATED TOOLS FOR SIGNAL RECOVERY

F E M T O