



Instrument Expert Original factory packaging www.dorgean.com

1.1.2.3 Low Power Thermal Sensors

1.1.2.3.1 Low Power BeamTrack-Power / Position / Size Sensors

100µW to 10W

Features (see introduction in pages 110-112)

- All the features of standard power sensors plus...
- Accurate tracking of beam position to fractions of a mm
- Monitoring of the laser beam size

3A-QUAD / 3A-P-QUAD





Model	3A-QUAD (a)	3A-P-QUAD (a)	10A-BB-16-PPS (a)
Use	General purpose	Short pulses	Low power
Functions	Power / Energy / Position	Power / Energy / Position	Power / Energy / Position / Size
Absorber Type	Low power broadband	P type	Broadband
Spectral Range µm	0.19 - 20	0.15 - 8	0.19 - 11
Aperture mm	Ø9.5mm	Ø12mm	Ø16mm
Power Mode			
Power Range	100µW - 3W	160µW - 3W	20mW - 10W
Power Scales	3W to 300µW	3W to 300µW	10W / 5W / 0.5W
Power Noise Level	5μW	10µW	1mW
Thermal Drift (30min)%	10 - 40µW ^(b)	10 - 40 μW ^(b)	NA
Maximum Average Power Density kW/cm ²	1	0.05	28
Response Time with Meter (0-95%) typ. s	1.8	2.5	0.8
Calibration Uncertainty ±%	1.9	1.9	1.9
Power Accuracy ±% ^(g)	3	3	3 ^(h)
Linearity with Power ±%	1	1	1
Energy Mode			
Energy Range	20µJ - 2J	30µJ - 2J	6mJ - 2J
Energy Scales	2J to 200µJ	2J to 200µJ	2J / 200mJ
Minimum Energy	20µJ	30µJ	6mJ
Maximum Energy Density J/cm ²			
<100ns	0.3	1 ^(f)	0.3
0.5ms	1	1 ^(f)	2
2ms	2	1 ^(f)	2
10ms	4	1 ^(f)	2
Beam Tracking Mode			
Position			
Beam Position Accuracy mm (c)	0.15	0.15	0.15
Beam Position Resolution mm	0.02	0.02	0.02
Min Power for Position Measurement	300µW	400µW	50mW
Size ^(d)			
Size Accuracy (e)	NA	NA	\pm (5%+50µm) for centered beam
Size Range mm (4o beam diameter)	NA	NA	1.5 - 10
Min Power for Size Measurement	NA	NA	50mW
Cooling	Convection	Convection	Convection
Weight kg	0.3	0.3	0.3
Fiber Adapter Available (see page 119)	ST, FC, SMA, SC	ST, FC, SMA, SC	ST, FC, SMA, SC
Compliance	CE, UKCA, China RoHS	CE, UKCA, China RoHS	CE, UKCA, China RoHS
Version	,,,		V1
Part number	7Z07934	7Z07935	7Z07905
Note: (a) The BeamTrack features are supported by Centauri, StarBright, StarLite, Nova II and Vega meters, Juno, Juno+, Juno-RS and EA-1 interfaces and StarLab application. Position and Size measurements work only in Power mode (but not in single shot Energy mode).			
Note: (b) Depending on room airflow and temperature val Note: (c) For position within inner 30% of aperture. Positi other desired position with Centauri, StarBright or StarLa	on measuring center corresponds to ge		ter can be software reset to geometric center or

Note: (d) Assumes laser beam with circular Gaussian (TEM₀₀) distribution. For other modes, size measurement is relative. Note: (e) Accuracy spec will be maintained for beams ≥1.8 mm not deviating from center by more than 15% of beam diameter. Note: (f) For P type and shorter wavelengths derate maximum energy density as follows: Wavelength Derate to value 1064nm not derated

vavelengtri	Derate to value
1064nm	not derated
532nm	not derated
355nm	40% of stated value
266nm	10% of stated value
193nm	10% of stated value
rated spectra	al curve at all wavelength

Note: (g) The 3A-QUAD has a relatively large spectral variation in absorption and has a calibrated spectral curve at all wavelengths in its spectral range to the above specified accuracy. Nova, Orion and LaserStar meters do not support this feature and when used with those meters, the accuracy will be ±3% as above for 532nm, 905nm, 1064nm and 10.6µm but there will be an additional error of up to 3% at other wavelengths in the spectral range 190 – 3000nm. Note: (h) ±4%. For wavelengths <240nm

* For drawings please see page 56

For latest updates, please visit our website: www.ophiropt.com





7

Interface Module on cable

0.5M cable to smart connector 1.5M cable to sensor 0 0 43 0 0 8 £ ---

3A-QUAD / 3A-P-QUAD



10A-BB-16-PPS

